


ELECTRIC WHEELCHAIR CATALOG

FREEDOM TO MOVE ANYWHERE



Pulse Motion Pro 1



PRODUCT	SPECIFICATIONS
Model Name	Pulse Motion Pro 1
Colour	Silver
Net Weight	44 KG
Product Material	Iron
Battery	Lithium-ion battery 24V 20Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front/Rear wheel size	10/16 inches
Brake	Electromagnetic Brake

Pulse Motion Pro 2



PRODUCT	SPECIFICATIONS
Model Name	Pulse Motion Pro 2
Colour	Black
Net Weight	41 KG
Product Material	Iron
Battery	Lithium-ion battery 24V 20Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front/Rear wheel size	10/24 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Cruza



PRODUCT	SPECIFICATIONS
Model Name	Pulse Cruza
Colour	Sliver
Net Weight	47 KG
Product Material	Iron
Battery	Lithium-ion battery 24V 20Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	10/16 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Autofold Smartride 2



PRODUCT	SPECIFICATIONS
Model Name	Pulse Autofold Smartride 2
Colour	Carbon Fiber
Net Weight	26 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 10Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushless Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Smartride 1



PRODUCT	SPECIFICATIONS
Model Name	Pulse Smartride 1
Colour	Carbon Fiber
Net Weight	26 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 10Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Xtrion



PRODUCT	SPECIFICATIONS
Model Name	Pulse Xtrion
Colour	Carbon Fiber
Net Weight	30 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 6Ah*2
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Innovax



PRODUCT	SPECIFICATIONS
Model Name	Pulse Innovax
Colour	Carbon Fiber
Net Weight	29 KG
Product Material	Aluminum
Battery	Lithium-ion Battery 24V 12Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Aerodrive 1



PRODUCT	SPECIFICATIONS
Model Name	Pulse Aerodrive 1
Colour	Metallic red & Black
Net Weight	37 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 12Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Aerodrive 2



PRODUCT	SPECIFICATIONS
Model Name	Pulse Aerodrive 2
Colour	Metallic Blue & Black
Net Weight	36 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 10Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Aluminum Alloy)
Brake	Electromagnetic Brake

Pulse Joylite 1



PRODUCT	SPECIFICATIONS
Model Name	Pulse Joylite 1
Colour	Carbon Fiber
Net Weight	21 KG
Product Material	Aluminum
Battery	Lithium-ion battery 24V 10Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	7/12 inches(Magnesium Alloy)
Brake	Electromagnetic Brake

Pulse Joylite 2



PRODUCT	SPECIFICATIONS
Model Name	Pulse Joylite 2
Colour	Metallic Blue
Net Weight	21 KG
Product Material	Iron
Battery	Lithium-ion battery 24V 10Ah
Charging Time	Up to 6 hour
Maximum Speed / Motor	6km/h / Brushed Motor 250W*2
Drive Range	15-20km
The Front / Rear wheel size	8/12 inches(Plastic)
Brake	Electromagnetic Brake

*Engineering Healthcare.
Enabling Life.*

Zuvio Technologies Private Limited

1st Floor, Plot No. 646. 27th Main Road, HSR Layout,
Sec-1, Bengaluru 560102.

+91 9071 101 108 www.pulseio.in info@pulseio.in



**To know more
Scan QR Code**